

Store at +2 to +8°C

PRODUCT CHARACTERISTICS:

CHRONOLAB Calibrator is a human lyophilized serum. The concentrations and activities have been selected so as to ensure optimum calibration for use in manual and automatic analysers.

REAGENTS:

Human serum. Biological additives. Bacteriostatics agents. The concentrations / activities of the components are lot – specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS:

Components from human origin have been tested and found negative for the presence of HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

CALIBRATION:

Information available on request.

PREPARATION:

Reconstitute with 3.0 ml of distilled water. Mix thoroughly, avoiding foam forming. Bring to room temperature for about 30 min. before use. Improper handling and/or storage can affect results. Inaccurate reconstitution and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY:

The calibrator is stable until the expiration date on the label when stored tightly closed at 2 - 8 °C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial grown. Store calibrators tightly capped when not use.

After reconstitution is stable for :

- at 15 °C to 25 °C - 8 hours
- at 2 °C to 8 °C - 2 days
- at -25 °C to -15 °C - 2 weeks

Bilirubin (stored protected from light):

- at 15 °C to 25 °C - 8 hours
- at 2 °C to 8 °C - 1 days
- at -25 °C to -15 °C - 2 weeks

Acid Phosphatase:

- at 15 °C to 25 °C - 6 hours
- at 2 °C to 8 °C - 2 days
- at -25 °C to -15 °C - 2 weeks

PACKAGING

Ref: 101-0353

Cont.: 20 x 3 ml